

Wheelchair Footrests

Sitting

The highly reliable and adjustable design of CAPS II footrest is now available for fitting on a range of standard wheelchairs including the B10 and Uni. They are fitted quickly and easily by removing the existing footrest tubes and replacing with this special footrest assembly.

They offer a wide range of adjustments so it is possible to accommodate a wide range of foot positions. Adjustment is provided front to back, side to side, plantar and dorsi flexion, inversion and eversion. Available in a range of sizes, and with either leather or webbing footstraps.

Size Information

WHEELCHAIR FOOTREST SIZE CHART				
Footplate Size (Length x Width)	Front Footstrap Length (standard)	Rear Footstrap Length (standard)	Maximum Heel Width	
A 16 x 10.5	30	40	7.5	
B 20 x 10.5	40	50	7.5	
C 22.5 x 11.5	40	50	8.5	
D 25 x 12.5	50	60	9.5	

All dimensions are in cm.



YOUR PRESCRIPTION DETAILS

FOOTREST	Footplate Size	_____	OPTIONAL	Leather Footstraps	<input type="checkbox"/>
DETAILS	Front Footstrap Length	_____	PARTS	Heel Loop Outriggers	<input type="checkbox"/>
	Rear Footstrap Length	_____			

Please note: there may be a small charge for some items.

Please attach a copy of this sheet to your order.

If you have other special requirements, please contact us to discuss the options.

Compatibility

The wheelchair footrests can be used as a direct replacement for standard wheelchair footrests where the footrest drop tube diameter is 3/4". This includes a number of standard wheelchairs such as the B10 and Uni.

Additional Support

They can be used to provide a specific foot position where the support offered by the wheelchair itself is not sufficient.

Kit Contents

The Wheelchair Footrests comprise of the following items:

- Black anodised footplates with fluted rubber infill
- Replacement footrest drop tubes to suit wheelchair footrest hangers supporting 3/4" diameter footrest drop tubes.
- Webbing footstraps
- Heel Loops
- Footplate support tubing
- All fixings and fittings.

Other Accessories

Leather Footstraps are available as an option.

Heel Loop Outriggers are available to increase the Maximum Heel Width of the footplate by approximately 16mm per side.

Please contact us for up to date prices.

We can also supply a one piece footplates. Please contact us if you require further information.

Other Accessories

Depending on the footrest height and position, castor outriggers may be required to ensure adequate castor clearance.

As a company we are committed to promoting a 24 hour approach to postural management through education, design and manufacture. The design of our equipment and training is based on a thorough understanding of the theoretical and practical requirements of the management of posture for people with moderate to severe motor impairment.



CE This product is manufactured to comply with the Medical Devices Directive MDD93/42/EEC.



Active Design Ltd
68K Wyrley Road
Birmingham B6 7BN
Tel: (0121) 326 7506
Fax: (0121) 327 8807
Web: www.activedesign.co.uk