

WHEELCHAIR COMPATIBILITY CHART

CAPS II & MiniCAPS Seating Systems



	Type	MiniCAPS															CAPS II									
		A1	A2	A3	B1	B2	B3	C1	C2	C3	D1	D2	D3	E1	E2	E3	Nursery	Nursery+	Child	Junior	Junior+	Teenage	Teenage+	Adult	Adult+	
8/9 Variants, Uni, Access, Accent	M	-	-	-	13	13	13	15	15	15	15	15	15	15 ¹	15 ¹	15 ¹	13	13	15	15	15 ¹	17	17	17 ¹	17 ¹	
Action 3/4 (Invacare)	M	12	12	12	13	13	13	14	14	14	15	15	15	16	16	16	12	13	14	15	16	17	18	19	20	
B400/B500 (Otto Bock)																						✓	✓	✓	✓	
Chunc	T	USE CHUNC ADAPT BASE																								
Discovery (Otto Bock)	T	12	12	12	14	14	14	14	14	14	16	16	16	16	16	16	14	14	16	16	16	16	18	18	18 ¹	
Bugzi (Meru)	P	✓	✓	✓	✓	✓	✓	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
Hi-Lo Base (R82)	S	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	3	3	3	-	-	
Ibis XP (Handicare)	P									✓	✓	✓	✓	✓	✓			✓	✓	✓	✓	✓	✓	✓		
Jazz Easy S (Tendercare)	B	1	1	1	1	1	1	1	1	1	2	2	2	2	2	2	1	1	1	2	2	-	-	-	-	
Kimba Spring/Neo (Otto Bock)	B	1	1	1	1	1	1	1	1	1	2	2	2	2	2	2	1	1	1	2	2	-	-	-	-	
Mini Tilt (Tendercare)	T	12	12	12	12	12	12	12	12	12	14	14	14	-	-	-	12	12	12	14	-	-	-	-	-	
Mojo (SOS)	T	-	-	-	-	-	-	-	-	-	Midi	Midi	Midi	Midi	Midi	Midi	-	-	-	Midi	Midi	Maxi	Maxi	Maxi	Maxi	
Neo (SW Seating)	T	12	12	12	13	13	13	14	14	14	15	15	15	16	16	16	12	13	14	15	16	17	18	19	20	
Puma 20 (Handicare)	P	-	-	-	-	-	-	-	-	-	-	-	-	✓	✓	✓	-	-	-	-	✓	✓	✓	✓	✓	
Rea Azalea (Invacare)	T	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	44	44	44	44	49	49	
Salsa Zippie/M/R (Sunrise Medical)	P	12	12	12	12	12	12	14	14	14	15	15	15	15	15	15	14	14	14	16	16	16	17	18	20	
Snappi (Tendercare)	B	✓	✓	✓	✓	✓	✓	✓	✓	✓	2	2	2	2	2	2	✓	✓	2	2	-	-	-	-	-	
Sparky (Handicare)	P	✓	✓	✓	✓	✓	✓	✓	✓	✓	-	-	-	-	-	-	✓	✓	✓	-	-	-	-	-	-	
Spectra Blitz (Invacare)	P	12	12	12	14	14	14	14	14	14	-	-	-	-	-	-	12	12	14	-	-	-	-	-	-	
Spectra/XTR/Topan (Invacare)	P	-	-	-	-	-	-	15	15	15	15	15	15	15	15	15	-	-	15	15	15	17	17	17	17	
Tiger/Serval (R82)	B	1	1	1	1	1	1	1	1	1	1	1	1	-	-	-	1	1	1	-	-	-	-	-	-	
Tom Cross/Streetwise (Activate)	B	Std	Std	Std	Std	Std	Std	Std	Std	Std	Maxi	Maxi	Maxi	Maxi	Maxi	Maxi	Std	Std	Std	Maxi	Maxi					
Zippie TS ² (Sunrise Medical)	T	12	12	12	13	13	13	14	14	14	14	14	14	16	16	16	12	12	14	16	18	18	18	18	-	
Modulite Base- All models (Invacare)	P	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	15	15	17	17	19	19	21	

¹ Armrest outriggers are required. ² Only the rigid version is compatible. The folding version can be used but the interface is rigidly attached and will prevent folding. Other bases are available - please contact us for further details. Base types: Buggy (B), Manual (M), Powered (P), Tilting (T)